

## REPLACEMENT PARTS

- Batteries
- Ear Cushions for Nardalerts
- Instrument Cases

### BATTERIES FOR NARDA METERS, SMARTS™, NARDALERT™

Instrument Model Number	Description	Part/Model Numbers		
		Narda	NEDA	IEC
SMARTS™ Circuit Fail-Safe	9.0V Lithium <sup>a</sup> 3.0V Lithium	88020000 88025000	1604LC —	— Note b
Nardalert™ Circuit LED <sup>d</sup>	12.0V Alkaline 1.5V Alkaline	88021000 88022000	1811A 135SO	Note c SR41
8699	6.0V Carbon	88016000	908AC	—
8217, 8512, 8520, 8711	9.0V Alkaline	88017000	1604A	6LR61
8616	+12.5V, -12.5V NiCad <sup>g</sup>	88011000	—	—
8716 <sup>d</sup> , 8719 <sup>d</sup> , 8850, 8870	8.4V NiCad <sup>g</sup>	88018000	—	—
8718, 8718B	8.4V NiCad <sup>g</sup>	21760000	—	—
8745, 8746, 8747, 8748	15.0V NiCad <sup>g</sup>	88024000	—	Note e
8211 <sup>d,f</sup> , 8611 <sup>d,f</sup>	6.0V Alkaline	88012000	1410M	4MR52
8110B	+15V, -15V NiCad <sup>g</sup>	88010000	—	—

#### NOTES:

- <sup>a</sup> 9.0V Alkaline can be substituted but must be replaced more often.  
<sup>b</sup> Renata 320A  
<sup>c</sup> Duracell MN21, Eveready A23, Ronson VR22, and GP23A  
<sup>d</sup> Requires two batteries  
<sup>e</sup> Motorola NLN-4462B  
<sup>f</sup> Duracell PC164A  
<sup>g</sup> Dispose of properly

### REPLACEMENT EAR CUSHIONS FOR NARDALERT™ PERSONAL MONITORS

The nardalert series of personal monitors are supplied with an earphone assembly that features quickly replaceable ear cushions. This permits multiple personnel to share the personal monitors yet still maintain hygiene standards. These ear cushions are sold twenty four pieces to a box.

Applicable Model Numbers	Narda Part Number*
All 884XB, 884XC, 884XD, and 884XE Series Models	86104003

\*Note: The part number for an individual ear cushion is 86104000 as shown in the nardalert manual. Part number 86104003 refers to a box of 24 pieces.

### INSTRUMENT CASES

Meter Model	Narda Part Number
8718, 8718B	32542707
8712, 8715	32542709
8511, 8513	32542708
8716, 8719, 8616, 8711	32542704

Meter Model	Narda Part Number
8512, 8520	32542702
8850, 8850B	32542701
8870	32542700
8110B	32542705